


Product Datasheet

PIP antibody (orb19224)

Description	Goat polyclonal antibody to PIP
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	GCDFP15 / PIP
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.IHC: In paraffin embedded Human Breast shows staining of luminal epithelial cells of the lactiferous duct. Recommended concentration, 2-3ug/ml.WB: Approx 14kDa band observed in lysates of cell line A431 (calculated MW of 16.6kDa according to NP_002643.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	GANKAQDNTRK
Clonality	Polyclonal
MW	16.6
Entrez	5304
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml, IHC-P: 2-3 µg/ml
Expiration Date	12 months from date of receipt.



orb19224
(0.3ug/ml)
staining of
A431 cel...



orb19224
(2.5ug/ml)
staining of
paraffin...