



Product Datasheet

PGRMC2 antibody (orb19227)

Description	Goat polyclonal antibody to PGRMC2
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF
Target	PGRMC2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:1000.WB: Preliminary experiments in human ovary, testes and placenta lysates gave no specific signal but low background (at antibody concentration up to 1 µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	DEEDTKDHNKQD
Clonality	Polyclonal
MW	26.2
Entrez	10424
Dilution Range	ELISA: 1:1000, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of paraforma...



Flow cytometric analysis of paraformalde...