

Product Datasheet

Carboxypeptidase A1/A2/B Rabbit mAb Antibody (orb1922724)

Description	Carboxypeptidase A1/A2/B Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	A synthesized peptide derived from human Carboxypeptidase A1+A2+B
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	Calculated MW: 47 kDa; Observed MW: 44 kDa
Uniprot ID	P15085, P48052, P15086
Dilution Range	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200
Expiration Date	12 months from date of receipt.