

Product Datasheet

Pyruvate Dehydrogenase E2 Rabbit mAb Antibody (orb1922727)

Description	Pyruvate Dehydrogenase E2 Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IP, WB
Immunogen	A synthesized peptide derived from human DLAT
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	Calculated MW: 69 kDa; Observed MW: 69 kDa
Uniprot ID	P10515
Dilution Range	WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/50 FC: 1/50-1/100
Expiration Date	12 months from date of receipt.