



Product Datasheet

PERP antibody (orb19231)

Description	Goat polyclonal antibody to PERP
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	PERP
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 20kDa band observed in human liver lysates (calculated MW of 21.4kDa according to NP_071404.2). Recommended concentration: 1-3 µg/ml.
Protein Sequence	CLPNYEDDLLGNAK
Clonality	Polyclonal
MW	21.4
Entrez	64065
Dilution Range	ELISA: 1:64000, WB: 1-3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human liver lys...



Immunohistochemical staining of Human Sk...