



Product Datasheet

Alpha E-catenin antibody (orb19232)


Description	Goat polyclonal antibody to Alpha E-catenin antibody
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	Catenin alpha1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 100kDa band observed in human colon and human skin lysates (calculated MW of 100kDa according to NP_001894.2). Recommended concentration: 1-3 µg/ml.
Protein Sequence	KREKQDETQTKIK
Clonality	Polyclonal
MW	100
Entrez	1495
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human colon lys...



Western blot analysis of HEK293 lysate u...



Immunohistochemical staining of Human Pr...