

## Product Datasheet

### CHD3 Rabbit mAb Antibody (orb1923262)

|                            |                                                                                                    |
|----------------------------|----------------------------------------------------------------------------------------------------|
| <b>Description</b>         | CHD3 Rabbit mAb Antibody                                                                           |
| <b>Species/Host</b>        | Rabbit                                                                                             |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                  |
| <b>Conjugation</b>         | Unconjugated                                                                                       |
| <b>Tested Applications</b> | FC, ICC, IHC-P, WB                                                                                 |
| <b>Immunogen</b>           | A synthesized peptide derived from human CHD3                                                      |
| <b>Preservatives</b>       | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| <b>Form/Appearance</b>     | Liquid                                                                                             |
| <b>Storage</b>             | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.           |
| <b>Note</b>                | For research use only                                                                              |
| <b>Isotype</b>             | IgG                                                                                                |
| <b>Clonality</b>           | Monoclonal                                                                                         |
| <b>MW</b>                  | Calculated MW: 227 kDa; Observed MW: 260 kDa                                                       |
| <b>Uniprot ID</b>          | <b>Q12873</b>                                                                                      |
| <b>Dilution Range</b>      | WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 FC: 1/50-1/100                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                    |