



Product Datasheet

AIRE antibody (orb19236)

Description	Goat polyclonal antibody to AIRE
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	AIRE
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.IHC: In paraffin embedded Human Thymus shows nuclear staining in select thymocytes. Recommended concentration, 2.5ug/ml.WB: Appr. 60kDa in transfected HEK293 lysates (calculated MW of 57.7kDa according to NP_000374.1). Recommended concentration: 0.5-1.0 µg/ml.
Protein Sequence	KDVDLSQPRKGRKP
Clonality	Polyclonal
MW	57.7; 36.5
Entrez	326
Dilution Range	ELISA: 1:32000, WB: 0.5-1.0 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Western blot analysis of HEK293 cell lys...