



Product Datasheet

AKR1C4 antibody (orb19237)

Description	Goat polyclonal to AKR1C4
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	AKR1C4
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody may cross-react with AKR1C1ELISA: Peptide ELISA: antibody detection limit dilution 1:4000.IHC: In paraffin embedded Human Liver shows staining of distinct structures in the cytoplasm of hepatocytes. Recommended concentration, 2-4ug/ml.WB: Approx. 35kDa band observed in Human Liver lysates (calculated MW of 37.1kDa according to NP_001809.2). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	DPKYQRVELNDGH
Clonality	Polyclonal
MW	37.1
Entrez	1109
Dilution Range	ELISA: 1:4000, WB: 0.1-0.3 µg/ml, IHC-P: 2-4µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human liver lys...



Immunohistochemical staining of Human Li...