

Product Datasheet

[KO] CD147 Mouse mAb Antibody (orb1924214)

Description	[KO] CD147 Mouse mAb Antibody
Species/Host	Mouse
Reactivity	Human, Monkey, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	Full length human recombinant protein of human BSG produced in 293T cell.
Preservatives	Liquid in PBS containing 50% glycerol, 1% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	0.64 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	Calculated MW: 42 kDa; Observed MW: 55 kDa
Uniprot ID	P35613
Dilution Range	WB: 1/500-1/1000 IHC: 1/50-1/100 ICC: 1/100-1/200
Expiration Date	12 months from date of receipt.