

Product Datasheet

[KO] STAT1 Mouse mAb Antibody (orb1924216)

Description	[KO] STAT1 Mouse mAb Antibody
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Full length human recombinant protein of human STAT1 produced in 293T cell.
Preservatives	Liquid in PBS containing 50% glycerol, 1% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	0.7 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	Calculated MW: 87 kDa; Observed MW: 87 kDa
Uniprot ID	P42224
Dilution Range	WB: 1/500-1/1000 ICC: 1/100-1/200
Expiration Date	12 months from date of receipt.