

Product Datasheet

[KO] SOD2 Mouse mAb Antibody (orb1924220)

Description	[KO] SOD2 Mouse mAb Antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human SOD2 produced in 293T cell.
Preservatives	Liquid in PBS containing 50% glycerol, 1% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	0.65 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	Calculated MW: 25 kDa; Observed MW: 20 kDa
Uniprot ID	P04179
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.