

Product Datasheet

[KO] CHCHD3 Mouse mAb Antibody (orb1924247)

Description	[KO] CHCHD3 Mouse mAb Antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Purified recombinant human CHCHD3 protein fragments expressed in E.coli.
Preservatives	Liquid in PBS containing 50% glycerol, 1% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	Calculated MW: 26 kDa; Observed MW: 26 kDa
Uniprot ID	Q9NX63
Dilution Range	WB: 1/500-1/1000 IHC: 1/50-1/100
Expiration Date	12 months from date of receipt.