

Product Datasheet

[KO] PUS7 Mouse mAb Antibody (orb1924248)

Description	[KO] PUS7 Mouse mAb Antibody
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human PUS7 produced in 293T cell.
Preservatives	Liquid in PBS containing 50% glycerol, 1% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	0.6 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	Calculated MW: 75 kDa; Observed MW: 90 kDa
Uniprot ID	Q96PZ0
Dilution Range	WB: 1/500-1/1000 FC: 1/50-1/100
Expiration Date	12 months from date of receipt.