

Product Datasheet

CD31 (5E2) Mouse mAb Antibody (orb1924455)

Description	CD31 (5E2) Mouse mAb Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Purified recombinant human CD31 protein fragments expressed in E.coli.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	4.5 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P16284
Dilution Range	IHC: 1/50-1/100
Expiration Date	12 months from date of receipt.