


Product Datasheet

MSN antibody (orb19245)

Description	Goat polyclonal antibody to MSN
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Moesin
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 70kDa band observed in Human Spleen lysates (calculated MW of 67.8kDa according to NP_002435.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	NARDESKKTAND
Clonality	Polyclonal
MW	67.8
Entrez	4478
Dilution Range	ELISA: 1:16000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb19245
(1µg/ml)
staining of
Human Sple...