

Product Datasheet

Estrogen Receptor alpha Mouse mAb Antibody (orb1924557)

Description	Estrogen Receptor alpha Mouse mAb Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Synthesized peptide derived from human Estrogen Receptor α
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P03372
Dilution Range	IHC: 1/100-1/500
Expiration Date	12 months from date of receipt.