

## Product Datasheet

### Cytokeratin 17 Mouse mAb Antibody (orb1924561)

|                            |                                                                                          |
|----------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>         | Cytokeratin 17 Mouse mAb Antibody                                                        |
| <b>Species/Host</b>        | Mouse                                                                                    |
| <b>Reactivity</b>          | Human                                                                                    |
| <b>Conjugation</b>         | Unconjugated                                                                             |
| <b>Tested Applications</b> | IHC-P                                                                                    |
| <b>Immunogen</b>           | Synthesized peptide derived from human Cytokeratin 17                                    |
| <b>Preservatives</b>       | Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.     |
| <b>Form/Appearance</b>     | Liquid                                                                                   |
| <b>Concentration</b>       | 1 mg/ml                                                                                  |
| <b>Storage</b>             | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                                    |
| <b>Isotype</b>             | IgG1                                                                                     |
| <b>Clonality</b>           | Monoclonal                                                                               |
| <b>Uniprot ID</b>          | <a href="#">Q04695</a>                                                                   |
| <b>Dilution Range</b>      | IHC: 1/100-1/500                                                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                          |