


Product Datasheet

PRNP antibody (orb19246)

Description	Goat polyclonal antibody to PRNP
Species/Host	Goat
Reactivity	Bovine, Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Prion Protein (143-153)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 25kDa band observed in Human Brain lysates (calculated MW of 27.7kDa according to NP_000302.1 and NP_898902.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	SDYEDRYREN
Clonality	Polyclonal
MW	27.7
Entrez	5621
Dilution Range	ELISA: 1:64000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



orb19246
(0.3ug/ml)
staining of
Human Br...