

Product Datasheet

GFP Tag Rabbit pAb Antibody (orb1924934)

Description	GFP Tag Rabbit pAb Antibody
Species/Host	Rabbit
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	A synthesized peptide derived from green fluorescent protein
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	Calculated MW: 27 kDa; Observed MW:
Uniprot ID	Q963I9
Dilution Range	WB: 1/1000-1/10000 IHC: 1/50-1/200 IF: 1/100-1/500
Expiration Date	12 months from date of receipt.