



Product Datasheet

MAP3K7 antibody (orb19250)

Description	Goat polyclonal antibody to MAP3K7
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, WB
Target	TAK1 / MAP3K7
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:1000.WB: Approx 75kDa band observed in lysates of cell lines NIH/3T3 and A431 (calculated MW of 64.2kDa according to NP_003179.1 and 67.2 according to NP_663304.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	AELDQDEKDQQNT
Clonality	Polyclonal
MW	64.2; 67.2
Entrez	6885
Dilution Range	ELISA: 1:8000, WB: 1-3 µg/ml, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of paraformaldehyde...



Western blot analysis of HeLa (lane 1)

a...