

Product Datasheet

Hepatitis C virus (HCV) core/NS3/NS4/NS5 protein (capture) (orb1925175)

Description	This recombinant protein contains immunodominant regions of HCV core, NS3, NS4, and NS5 proteins wit...
Tested Applications	CLIA, ELISA, LF, WB
Preservatives	0.1% ProClin300. Phosphate buffer saline, pH 8.0, and 8 M Urea
Form/Appearance	Purified, Liquid
Storage	Upon receipt, store at 2-8°C
Note	For research use only
Purity	> 90% pure by SDS-PAGE
Expiration Date	6 months from date of receipt.