

Product Datasheet

Hepatitis C virus (HCV) NS3 protein (AA1192-1457) (orb1925178)

Description	This recombinant protein contains immunodominant regions of HCV NS3 protein with a MW of 30.3 kDa. I...
Tested Applications	CLIA, ELISA, LF, WB
Preservatives	0.1% ProClin300. Phosphate buffer saline, pH 8.0, and 8 M Urea
Form/Appearance	Purified, Liquid
Storage	Upon receipt, store at 2-8°C
Note	For research use only
Purity	> 90% pure by SDS-PAGE
Expiration Date	6 months from date of receipt.