

Product Datasheet

Hepatitis C virus (HCV) NS3 protein (AA1192-1367) (orb1925180)

Description This recombinant protein contains immunodominant regions of HCV NS3 protein with a MW of 44.5 kDa. I...

Tested Applications CLIA, ELISA, LF, WB

Preservatives 0.1% ProClin300. Tris-buffered saline, pH 8.0

Form/Appearance Purified, Liquid

Storage Upon receipt, quick spin before open vial, aliquot and store at $\leq -20^{\circ}\text{C}$. Avoid multiple freeze/thaw cycles

Note For research use only

Purity > 90% pure by SDS-PAGE

Expiration Date 6 months from date of receipt.