

Product Datasheet

SCAP Antibody (N-term) (orb1925366)

Description	Purified Mouse Monoclonal Antibody (Mab)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This SCAP antibody is generated from a mouse immunized with a KLH conjugated synthetic peptide between 49-83 amino acids from the N-terminal region of human SCAP.
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1, κ
Clonality	Monoclonal
MW	139729 Da
Uniprot ID	Q12770
Dilution Range	WB: 1:4000
Expiration Date	12 months from date of receipt.