

Product Datasheet

TAS2R50 Antibody (C-Term) (orb1925657)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This TAS2R50 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 198-222 amino acids from human TAS2R50.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	34558 Da
Uniprot ID	P59544
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.