



Product Datasheet

CHGA antibody (orb19258)

Description	Goat polyclonal antibody to CHGA
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, WB
Target	Chromogranin A precursor
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:4000.WB: Approx 50kDa band observed in human kidney lysates (calculated MW of 50.7 kDa according to NP_001266.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	QGAKERAHQQKKH
Clonality	Polyclonal
MW	50.7; 34.3
Entrez	1113
Dilution Range	ELISA: 1:4000, WB: 1-2 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



orb19258 (0.1µg/ml)
staining of human
ki...



Immunohistochemical
staining of Human
Sm...