



Product Datasheet

IGSF9 antibody (orb1926)


Description	IGSF9 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human IGSF9 (551-680/1179 aa)
Target	IGSF9
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	124 kDa
Uniprot ID	Q9P2J2
Dilution Range	Flow-Cyt=1µg /Test
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of Hepg2 Cells ...



Immunohistochemical staining of rat kidn...



Immunohistochemical staining of human he...