

Product Datasheet

PLAGL1 Antibody (N-Term) (orb1926002)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This PLAGL1 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 112-146 amino acids from human PLAGL1.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	50819 Da
Uniprot ID	Q9UM63
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.