

Product Datasheet

UQCRFS1P1 Antibody (Center) (orb1926037)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This UQCRFS1P1 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 71-104 amino acids from the Central region of human UQCRFS1P1.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	30815 Da
Uniprot ID	P0C7P4
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.