

Product Datasheet

Mouse Noto Antibody (C-term) (orb1926625)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	This Mouse Noto antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 200-234 amino acids from the C-terminal region of Mouse Noto.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	26322 Da
Uniprot ID	Q5TIS6
Dilution Range	WB: 1:2000, WB: 1:4000
Expiration Date	12 months from date of receipt.