

Product Datasheet

NFIA Antibody (C-term) (orb1926778)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Target	This NFIA antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 384-418 amino acids from the C-terminal region of human NFIA.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	55944 Da
Uniprot ID	Q12857
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.