

Product Datasheet

(Mouse) Sox2 Antibody (N-term) (orb1926835)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	This Mouse Sox2 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 6~39 amino acids from the N-terminal region of mouse sox2.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	34454 Da
Uniprot ID	P48432
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.