

Product Datasheet

CILP antibody (orb19269)

Description	Goat polyclonal antibody to CILP
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CILP
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000. WB: Approx 150kDa band observed in rat spinal cord lysates (calculated MW of 133kDa according to human NP_003604.2 and of 132kDa according to Rat NP_001101631.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	NYRRTDHEDPRVK
Clonality	Polyclonal
MW	133
Entrez	8483
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



orb19269
(0.03ug/ml)
staining of rat
spi...