

Product Datasheet

Myoglobin Antibody (C-term) (orb1927055)

Description	Purified Mouse Monoclonal Antibody (Mab)
Species/Host	Mouse
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	This Myoglobin antibody is generated from a mouse immunized with a KLH conjugated synthetic peptide between 118-150 amino acids from the C-terminal region of human Myoglobin.
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	17184 Da
Uniprot ID	P02144
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.