

Product Datasheet

M Tcfap2d Antibody (C-term) (orb1927328)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	This Mouse Tcfap2d antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 308-335 amino acids from the C-terminal region of mouse Tcfap2d.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	49548 Da
Uniprot ID	Q91ZK0
NCBI	NP_694794.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.