


Product Datasheet

POU5F1 antibody (orb19274)

Description	Goat polyclonal antibody to POU5F1
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	OCT4 / POU5F1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 38kDa band observed in Human Liver lysates (calculated MW of 38.6kDa according to NP_002692.2). Recommended concentration: 1-3 µg/ml.
Protein Sequence	C-VTTLGSPMHSN
Clonality	Polyclonal
MW	38.6; 30.1
Entrez	5460
Dilution Range	ELISA: 1:2000, WB: 0.5--2 µg/ml
Expiration Date	12 months from date of receipt.



orb19274
(2µg/ml)
staining of
Human Pla...



orb19274
(0.5µg/ml)
staining of
Mouse H...