


Product Datasheet

TNS1 antibody (orb19275)

Description	Goat polyclonal antibody to TNS1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	Tensin 1 / TNS1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.WB: Preliminary experiments in Human Heart, Kidney and Skeletal Muscle lysates gave no specific signal but low background (at antibody concentration up to 1 µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	VSKVMLNAGQKR
Clonality	Polyclonal
MW	186
Entrez	7145
Dilution Range	ELISA: 1:2000, WB: 1 µg/ml, IHC-P: 4-6µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Kidney tissue showing TNS1 expression.