



Product Datasheet

ROBO1 antibody (orb19276)

Description	Goat polyclonal antibody to ROBO1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Target	ROBO1 / DUTT1 (Internal)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 250kDa band observed in Human Endothelial Cell line HUVEC lysates (calculated MW of 181kDa according to NP_002932.1). Some of the observed proteolytic breakdown products are also found in Human, Mouse and Rat Brain lysates. Recommended concentration: 1-3 µg/ml.
Protein Sequence	GDVDLSNKINEMK
Clonality	Polyclonal
MW	181; 177
Entrez	6091
Dilution Range	ELISA: 1:8000, WB: 1 µg/ml, IHC-P: 2 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HUVEC lysate us...



Immunohistochemical staining of Human Li...