

Product Datasheet

P2R3C Antibody (N-term) (orb1927754)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	This P2R3C antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 50-79 amino acids from the N-terminal region of human P2R3C.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	53316 Da
Uniprot ID	Q969Q6
NCBI	NP_060387.2
Dilution Range	WB: 1:1000, IHC-P: 1:50~100
Expiration Date	12 months from date of receipt.