

Product Datasheet

Myosin IXB antibody (orb19278)

Description	Goat polyclonal to Myosin IXB
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Myosin IXB
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.WB: Approx 200kDa band observed in Human Liver lysates (calculated MW of 243kDa according to NP_004136.2). Recommended concentration: 2-6 µg/ml.
Protein Sequence	KNDFEPVKQSKAQ
Clonality	Polyclonal
MW	243; 229
Entrez	4650
Dilution Range	ELISA: 1:2000, WB: 2-6 µg/ml
Expiration Date	12 months from date of receipt.



orb19278
(2ug/ml)
staining of
Human
Live...