

## Product Datasheet

### HSP90B1 antibody (orb19279)

## Description

Goat polyclonal antibody to HSP90B1



## Species/Host

Goat

## Reactivity

Bovine, Canine, Human, Mouse, Porcine, Rat, Zebrafish

## Conjugation

Unconjugated

## Tested

ELISA, FC, IF, IHC, WB

## Applications

Western blot analysis of NIH3T3 cell lys...

## Target

GPR94 / TRA1

## Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

## Storage


Aliquot and store at -20°C. Minimize freezing and thawing.

## Note


For research use only

## Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:32000. IHC: In paraffin embedded Human Liver shows staining in some hepatocyte nuclei and weak cytoplasm staining of hepatocytes. Recommended concentration, 3-5ug/ml. WB: Approx 100kDa band observed in lysates of cell line NIH/3T3 (calculated MW of 92.5kDa according to NP\_003290.1). Recommended concentration: 0.01-0.03 µg/ml.



Immunohistochemical staining of Human Ut...



Western blot analysis of U2OS (lane 1), ...

## Protein Sequence

KEGVKFDESEKTKE

## Clonality

Polyclonal

## MW

92.5

## Entrez

**7184**

## Dilution Range

ELISA: 1:32000, WB: 0.1-0.3 µg/ml, IHC-P: 3-5 µg/ml

## Expiration Date

12 months from date of receipt.