

Product Datasheet

RCAS1 antibody (orb19280)

Description

Goat polyclonal antibody to RCAS1

Species/Host

Goat

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Target

EBAG9 / RCAS1

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:64000. WB: Approx 26kDa band observed in lysates of cell lines HEK293 and A431 (calculated MW of 24.4kDa according to NP_004206.1 and NP_004206.1). Recommended concentration: 1-3 µg/ml. A minor band of unknown identity was also consistently observed at 75kDa. This band was successfully blocked by incubation with the immunising peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

KEAQLMKKEQN

Clonality

Polyclonal

MW

24.4

Entrez

9166

Dilution Range


ELISA: 1:64000, WB: 1-3 µg/ml, IHC-P: 3.75 µg/ml

Expiration Date

12 months from date of receipt.



Western blot analysis of A431 lysate using...



Immunohistochemical staining of Human Sk...