

Product Datasheet

UPK1A Antibody (C-term) (orb1928032)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This UPK1A antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 191-217 amino acids from the C-terminal region of human UPK1A.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	28879 Da
Uniprot ID	O00322
NCBI	NP_001268372.1 , NP_008931.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.