


Product Datasheet

HTR1B antibody (orb19283)

Description	Goat polyclonal antibody to HTR1B
Species/Host	Goat
Reactivity	Bovine, Canine, Feline, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Serotonin receptor 1B / HTR1B
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 48kDa band observed in Human Brain lysates (calculated MW of 43.6kDa according to NP_000854.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	SNEDFKQAFHK
Clonality	Polyclonal
MW	43.6
Entrez	3351
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb19283
(2ug/ml)
staining of
Human Brai...



Western blot
analysis of
Human
Cerebellu...