


Product Datasheet

RNF115 antibody (orb19287)

Description	Goat polyclonal antibody to RNF115
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Rabring 7 / BCA2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 35kDa band observed in Human lysates of cell line HEPG2 (calculated MW of 33.7kDa according to NP_055270.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	DQDNRRANERGHQTH
Clonality	Polyclonal
MW	33.7
Entrez	27246
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HEPG2 cell lysa...