

Product Datasheet

MAP3K5 Antibody(T842) (orb1928771)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	This MAP3K5 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 821-849 amino acids from human MAP3K5.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	154537 Da
Uniprot ID	Q99683
NCBI	NP_005914.1
Dilution Range	WB: 1:2000, WB: 1:1000, IHC-P: 1:10~50
Expiration Date	12 months from date of receipt.