

## Product Datasheet

### MEKK4 Antibody (C-term) (orb1928774)

|                        |                                                                                                                                                                                           |
|------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Purified Rabbit Polyclonal Antibody (Pab)                                                                                                                                                 |
| <b>Species/Host</b>    | Rabbit                                                                                                                                                                                    |
| <b>Reactivity</b>      | Human                                                                                                                                                                                     |
| <b>Conjugation</b>     | Unconjugated                                                                                                                                                                              |
| <b>Target</b>          | This MEKK4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1074~1104 amino acids from the C-terminal region of human MEKK4.                  |
| <b>Form/Appearance</b> | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS. |
| <b>Note</b>            | For research use only                                                                                                                                                                     |
| <b>Isotype</b>         | Rabbit IgG                                                                                                                                                                                |
| <b>Clonality</b>       | Polyclonal                                                                                                                                                                                |
| <b>MW</b>              | 181685 Da                                                                                                                                                                                 |
| <b>Uniprot ID</b>      | <a href="#">Q9Y6R4</a>                                                                                                                                                                    |
| <b>NCBI</b>            | <a href="#">NP_005913.2</a> , <a href="#">NP_006715.2</a>                                                                                                                                 |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                                                                                                                           |