

Product Datasheet

MEK2 Antibody (N-term) (orb1928777)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Target	This MEK2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1-30 amino acids from the N-terminal region of human MEK2.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	44424 Da
Uniprot ID	P36507
NCBI	NP_109587.1
Dilution Range	IF: 1:10~50, WB: 1:1000, IHC-P: 1:50~100
Expiration Date	12 months from date of receipt.