

Product Datasheet

MEK1 (MAP2K1) Antibody (N-term) (orb1928778)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	This MEK1 (MAP2K1) antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1-30 amino acids from the N-terminal region of human MEK1 (MAP2K1).
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	43439 Da
Uniprot ID	Q02750
NCBI	NP_002746.1
Dilution Range	WB: 1:1000, IHC-P: 1:10~50
Expiration Date	12 months from date of receipt.